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P. 61 - P. 74

**ISO 9001:2008**

Ultrasonic Milk Analyzers EKOMILK are certified by the Bulgarian, Ukrainian, Czech, Rumanian, Polish, Russian and Belarusian Centers of Metrology.

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**BULTEH 2000 Ltd.**

Industrial Area, VEREA PLAST building

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e-mail: [bulteh@sz.inetg.bg](mailto:bulteh@sz.inetg.bg)

## ENGLISH

**EKOMILK** Ultrasonic Milk Analyzers are designed for fast and cost effective analysis of the fat contents, milk solids non fat (SNF), protein, milk density of cow, sheep, buffalo, goat, etc. milk and added water to milk as well.

### KEY FEATURE

- Simple and lightweight design;
- +12V DC and 220V (or 110V optional) AC Power supply;
- Cost effective:
  - Low power consumption;
  - Very small quantity of milk required;
  - No acid or other chemicals are used;
- One year full warranty;
- Measuring accuracy adjustment can be done by the user;
- RS 232 Interface;
- ESC POS Printer Support.

### ENVIRONMENTAL CONDITIONS:

Ambient air temperature	10° - 30°C
Milk temperature	10° - 30°C
Relative humidity	30% - 80%

### ELECTRICAL PARAMETERS:


AC Power Supply voltage	see p. 3
DC Power Supply voltage	12V to 14,2V
Power Consumption	30W max (100W - Ekomilk UltraPro)

### MECHANICAL PARAMETERS:

Dimensions (WxDxH)	95 x 300 x 250 mm
Weight	< 4.0 kg


### MEASURING PARAMETERS:

Fat	from 0,5% to 12% with accuracy $\pm 0,1\%$
Solids non fat (SNF)	from 6% to 12% with accuracy $\pm 0,2\%$
Milk density	from 1,0260 g/cm <sup>3</sup> to 1,0330 g/cm <sup>3</sup> $\pm 0,0005$ g/cm <sup>3</sup>
Protein	from 2% to 6% with accuracy $\pm 0,2\%$
Freezing point	from 0 to -1.000°C with accuracy $\pm 0,015^\circ\text{C}$
Added water to milk	from 0% to 60% with accuracy $\pm 5\%$

**CONTROLS  
(on the front panel)**


1. Select the work mode
2. Skip and search forwards and backwards
3. Confirm the choice
4. Tube (sucker)
5. Plastic support

**SWITCHES and OUTLETS  
(on the rear panel)**

1. RS 232 Interface outlet
2. 12V DC Printer Power supply outlet
3. ON/OFF Switch - for car vehicle use 12V
4. 12V DC Power supply outlet
5. POWER Switch
6. AC power led

## INSTRUCTIONS FOR USE

---

### ANALYZER INATALLATION

---

Place vertically the Analyzer on a table or any other flat surface.

- AC Power Supply Voltage

Connect the AC power lead to the mains socket.

Set the Power switch to "On".

- 12V DC Power Supply Voltage

Connect the DC supply lead to 12V outlet (black is "-") and to the autonomous DC supply (for example car supply).

Set the +12 V Power switch to "On".


#### Warm up

When the power is on, **WARM UP** appears on the display. When the "warm up" stage is over in about 5 minutes, **EKOMILK** is shown on the display. The Analyzer is ready to use.

Do not take into consideration the first sample because it is likely to be with a deviation out of specification. In case the analyzer has not been in use for more than a week (after transport for example), please first make 3 to 5 measurements of clean water without taking in consideration the measurement readings.

**Warning: The covers of the Analyzer should never be removed while the power leads are connected.**

**Warning: Under no circumstance should you try to repair the Analyzer's power lead yourself. In case of power lead damage, contact your dealer to make the repairs.**

---

### MILK SAMPLES

---

Milk samples temperature should be 10 - 30°C. If the milk temperature is above 38°C the message **HOT SAMPLE** appears on the display.

If you try to test cool (refrigerated) milk which has some milk fat/cream separated you will probably get wrong result especially for the milk fat contents. In this case you need to warm up the milk up to 40 - 42°C first, mix the milk in order to solve the separated fat, cool it down to 20 - 25°C and then you can test it with the EKOMILK Analyzer.

The acidity of the milk sample must be less than 25°T for cow, buffalo and goat milk and less than 28°T for sheep milk.

Use the milk sample only once. When the measuring is carried out, throw the sample away.

---

### BASIC MODES FOR USE

---

Press the **MODE** button to enter basic modes menu. Basic modes available are as follows:

- **COW MILK** - analysis of cow milk
- **CLEANING** \* - cleaning at the end of working day (Available in **MOTOR ON** mode)
- **MOTOR CHOICE\*** - choice the basic mode (Ekomilk, Ekomilk-120 and Ekomilk-M)
- **CALIBRATION** - calibration
- **SYSTEM** - manufacturer's mode only
- **SHEEP (BUFFALO, GOAT) MILK** - analysis of sheep milk ( buffalo or goat milk - see p. 3)

\* These modes are available for Milk analyzers with suction pump installed only.

Press the search buttons ▲, ▼ to select desired mode.

---

## MILK ANALYSIS WITH MOTOR (SUCTION PUMP)


---

Press the **MODE** button only once. Press the search buttons ▲, ▼ to select **MOTOR CHOICE** mode.

Press **OK** when the **MOTOR CHOICE** mode is displayed.

Set the desired motor mode - **ON** with search buttons ▲, ▼. Press **OK** to confirm it. The message **MOTOR TURNED ON** appears on the display. This way the “**MOTOR ON**” mode will remain set, even if the power supply is turned off.

Take out the plunger from the syringe (Step1) (for Ekomilk Ultra Pro - both plungers). Insert the plastic plug with the vinyl tube instead of the plunger (Step 2).

**The Analyzer is ready for use in mode with suction pump.**

### Measurement

**Steps:** Fill the measuring mug with the milk sample to be measured. Place the measuring mug on the plastic support (5) with the tube (sucker) into the milk sample.

Press **MODE** and by means of the search buttons ▲, ▼ select the desired mode:

- **COW MILK** - analysis of cow milk
- **SHEEP (BUFFALO, GOAT) MILK** - analysis of sheep milk (buffalo, goat etc. milk - see p. 3)

When the proper type of milk is displayed, press **OK** to start the measurement. The milk is automatically driven into the measuring camera. The message **WORKING** appears on the display while the measurement is started.

**Warning: Do not remove the measuring mug while the measurement is going on. This might cause result deviations out of the spec limits.**

**Warning: In case of formation of air bubble in the measuring camera the message EMPTY CAMERA will appear on the display. In this case make the measurement again.**

As soon as the measurement is completed the display shows the results for following milk parameters:

**FAT, SOLIDS NON FAT (SNF), MILK DENSITY** - on the upper row

**ADDED WATER TO MILK, FREEZING VALUE\*, PROTEIN** - on the bottom row

\* $FREEZING\ VALUE = FREEZING\ POINT\ [^{\circ}C] \cdot (-100)$

The milk returns automatically into the mug.

**Warning: In case of motor damage the message MOTOR ERROR will appear on the display. Contact your dealer to make the repairs. Meanwhile you may use the Analyzer in MOTOR OFF mode. (see "MILK ANALYSIS WITHOUT SUCTION PUMP")**

---

#### **PRINTING DATA (IN CASE OF MICRO PRINTER CONNECTED TO THE ANALYZER)**

---

After the current measurement is completed the results could be printed by pressing the arrow button **▲** on the analyzer front panel. Results are printed every time when this button is pressed.

---

#### **MILK ANALYSIS WITHOUT MOTOR (SUCTION PUMP) AVAILABLE FOR MILK ANALYZER TYPE EKOMILK, EKOMILK120, EKOMILK-M**

---

Press the **MODE** button once. Press the search buttons **▲**, **▼** to select the **MOTOR CHOICE** mode. Press **OK** when the **MOTOR CHOICE** mode is displayed. Set the desired motor mode - **OFF** with search buttons **▲**, **▼**. Press **OK** to confirm it. The message **MOTOR TURNED OFF** appears on the display. This way the "**MOTOR OFF**" mode will remain set, even if the power supply is turned off.

Take out the plastic plug with the vinyl tube from the syringe. Insert the plunger instead of the plastic plug with the vinyl tube.

**The Analyzer is ready for use in mode without suction pump.**

#### **Measurement**

**Steps:** Press **MODE** button once and by means of the search buttons **▲**, **▼** select the desired mode:

■ **COW MILK** - analysis of cow milk

■ **SHEEP (BUFFALO, GOAT) MILK** - analysis of sheep milk (buffalo, goat milk etc. - see p. 3)

When the proper type of milk is displayed, fill the measuring mug with milk sample to be measured. Place the measuring mug on the plastic support (5) with the tube (sucker) into the milk sample. Pull up very carefully the plunger. This way the milk is driven into the measuring camera.

**Warning:** *Pull up the plunger smoothly to prevent formation of air bubbles in the measuring camera. In case of formation of air bubbles in the measuring camera the message **EMPTY CAMERA** will appear on the display. In this case make the measurement again.*

Press **OK** to start the measurement.

The message **WORKING** appears on the display while measurement is going on.

**As soon as the measurement is completed the display shows the results for following milk parameters:**

**FAT, SOLIDS NON FAT (SNF), MILK DENSITY** - on the upper row

**ADDED WATER TO MILK, FREEZING VALUE\*, PROTEIN** - on the bottom row

\*FREEZING VALUE = FREEZING POINT [°C] . (-100)

Pull down the plunger to pump out the milk from the measuring camera into the mug and throw the used milk sample away.

---

## CLEANING OUT

---

**This section has for an object to give directions for daily and periodical cleaning of milk analyzers Ekomilk with a view to protect their measuring systems from contamination and ensure troublefree work for a long time.**

### 1. DAILY PROCEDURE

It is necessary to clean the analyzer within the working day in the following cases:

- The interval between two consecutive measurements is more than half an hour;
- End of working day.

#### 1.1. Cleaning out when the interval between two consecutive measurements is more than half an hour.

Steps:

1.1.1. Fill the measuring mug with clean, warm (but not hot) and clean water (40°- 60°C).

1.1.2. Press the **MODE** button once. Press the search buttons **▲**, **▼** to select **CLEANING** option. Confirm with **OK**. The display shows message **CYCLES 01** - number of cleaning cycles to be done. One cleaning cycle pumps the water in and out of the sensor five times. Press the search buttons **▲**, **▼** to set the desired number of cleaning cycles. This number can range from 1 to 99. We recommend one or two cleaning cycles to be chosen. When the desired number is displayed, press **OK** to start the cleaning. When the **CLEANING** stage is over **CLEANING END** is shown on the display.

1.1.3. Remove the mug and throw away the muddy water. Repeat this procedure several times till clean water comes out from the Analyzer.

1.1.4. Done.

### **1.2. Cleaning out at the end of a working day.**

This procedure prevents formation and collection of fat and “Milk stone” deposits into the sensor. Milk stone consists of milk solids, calcium, magnesium, iron, sulfates, etc. Milk and water mineral deposits become hardened and layered on the sensor and vinyl pipes inner surfaces, which contact with milk and disturbs the milk analyzer work. Cleaning will be effective if a reagent that attacks the “milk stone” is used. We recommend 2% solution of the alkaline cleaner **EkoDay** to be used as a daily cleaning solution.

Preparation of 2% **EkoDay** working solution:

- a. Use pipette to add 10 mL **EkoDay** to glassware with 490 mL distilled water.
- b. Pour into a labeled container.

Take care this solution does not contact your eyes or skin!

**Attention! Use only 2% EkoDay working solution as a cleaning agent. The EkoDay concentrate can damage your analyzer!**

Steps:

1.2.1. First clean the analyzer with pure water as it is described in procedure 1.1. (see above).

**Attention! Using 2% solution of the alkaline cleaner EkoDay without first cleaning analyzer for removing the fats and proteins will result in fixing the “milk stone” to the surface. Always clean analyzer with pure water before using the cleaning agent.**

1.2.2. Fill the measuring mug with 2 % solution of the alkaline cleaner **EkoDay** (25°- 40°C), put the analyzer in cleaning mode, set 20 cycles and press the **OK** button. When the **CLEANING** stage is over **CLEANING END** is shown on the display.

1.2.3. Fill the measuring mug with clean water, put the analyzer in cleaning mode, set 5 cycles and press the **OK** button. When this procedure is done, remove the mug and throw away the water. Fill the measuring mug with clean water and repeat this procedure 3-4 times.

1.2.4. Done

## **2. PERIODICAL PROCEDURE**

To ensure a good work of the Milk Analyzer it is advisory to clean the device at least once a week strictly performing underwritten procedure.

This procedure uses 5% solution of acid cleaner **EkoWeek** as a periodical cleaning solution.

Preparation of 5% **EkoWeek** working solution:

- a. Use pipette to add 25 mL **EkoWeek** to glassware with 475 mL distilled water.
  - b. Pour into a labeled container.
- Take care this solution does not contact your eyes or skin!

**Attention! Use only 5 % EkoWeek working solution as a cleaning agent. The EkoWeek concentrate can damage your analyzer!**

2.1. Clean the analyzer first as it is cleaned at the end of the working day;

**Attention! Make sure that the analyzer is properly cleaned with pure water before going on to the next item of the procedure. Mixing both cleaning solutions- alkaline EkoDay and acid EkoWeek will result in forming “milk stone”.**

2.2. Fill the measuring mug with 5 % solution of the acid cleaner **EkoWeek** (25°- 40°C), put the analyzer in cleaning mode, set 40 cycles and press the **OK** button. When the **CLEANING** stage is over, **CLEANING END** is shown on the display.

2.3. Take the plastic plug (for Ekomilk UltraPro both plastic plugs) with the vinyl tube out of the holes. Insert the plunger instead of the plastic plug with the vinyl tube.

2.4. Fill the measuring mug with clean and warm, but not hot, water (40°- 60°C).

Pull up and down the plunger several times. Remove the mug and throw away the water. Fill the mug with clean and warm water and repeat this step 4-5 times.

2.5. Take the plunger out of the syringe. Wait until all water comes out. Insert back slowly the plunger to the bottom.

2.6. Remove the measuring mug. Take the plunger out of the syringe again. This time insert the rubber plug with the vinyl tube instead of the plunger.

2.7. Done

**Attention!**

1. Perform strictly this procedure to remove both acid cleaning solution and milk stone from the analyzer measuring system.
2. Do not use hard water (water that contains high amounts of  $\text{Ca}^{2+}$  and/ or  $\text{Mg}^{2+}$ ) for analyzer cleaning. For best results use distilled or deionized water.

---

## CALIBRATION

---

The Analyzer should be calibrated if the measuring accuracy for one or more milk parameters is out of the specified limits.

### **Calibration values determination:**

**- Determine the milk sample FAT contents by means of a classical method (for example you can use Gerber method);**

- Test the same milk using Milk analyzer;

- Subtract the Milk analyzer **FAT** value from **FAT** value achieved with classical method. If the difference value is in the spec limits there is no need of fat calibration. Otherwise add this difference value as a fat calibration value using the procedure described in "**Calibration values Saving**";

**- Determine the milk sample SOLIDS NON FAT (SNF) by means of a classical method;**

- Test the same milk using Milk analyzer;

- Subtract the Milk analyzer **SOLIDS NON FAT (SNF)** value from **SOLIDS NON FAT (SNF)** value achieved with classical method. If the difference value is in the spec limits there is no need of **SOLIDS NON FAT (SNF)** calibration. Otherwise add this difference value as a **SOLIDS NON FAT (SNF)** calibration value using the procedure described in "**Calibration values saving**";

**- Determine the milk sample DENSITY by means of a density meter;**

- Test the same milk using Milk analyzer;

- Subtract the Milk analyzer **DENSITY** value from **DENSITY** value achieved with density meter. If the difference value is in the spec limits there is no need of milk **DENSITY** calibration. Otherwise add this difference value as a milk **DENSITY** calibration value using the procedure described in "**Calibration values Saving**";

**Note: If you will calibrate both SOLIDS NON FAT (SNF) and DENSITY parameters, first calibrate SOLIDS NON FAT (SNF). Test the same milk using EKOMILK Analyzer again. Calibrate the DENSITY only if it is necessary.**

**- Determine the milk sample PROTEIN contents by means of a classical method (for example you can use Kjeldahl method);**

- Test the same milk using EKOMILK Analyzer;

- Subtract the Milk analyzer **PROTEIN** value from **PROTEIN** value achieved with classical method. If the difference value is in the spec limits there is no need of protein calibration. Otherwise add this difference value as a protein calibration value using the procedure described in "**Calibration values Saving**";

**- Determine the milk sample FREEZING POINT by means of a classical method (for example you can use Cryoscope);**

- Test the same milk using EKOMILK Analyzer;

- Subtract the Milk analyzer **FREEZING POINT** value from **FREEZING POINT** value achieved with classical method. If the difference value is in the spec limits there is no need of freezing point calibration. Otherwise add this difference value as a freezing point calibration value using the procedure described in "**Calibration values Saving**";

*Example 1:*

*Freezing point by classical method - (-0.548)*

Freezing value by EKOMILK - (53.0)  
 Freezing point by EKOMILK - (-0.53)  
 Freezing value correction value= $(-0.548) - (-0.53) = (-0,018)$   
 Example2:  
 Freezing point by classical method - (-0.548)  
 Freezing value by EKOMILK - (56.0)  
 Freezing point by EKOMILK - (-0.56)  
 Freezing value correction value= $(-0.548) - (-0.56) = (+0,012)$

**Attention:**

1. If you need to calibrate both **SOLIDS NON FAT (SNF)** and **FREEZING POINT**, first calibrate **SOLIDS NON FAT (SNF)**. Test the same milk using **EKOMILK Analyzer** again. Calibrate the **FREEZING POINT** only if it is necessary.
2. The freezing point correction will change also the Added water measuring result.

---

## CALIBRATION VALUES SAVING

---

- To save a new calibration value press **MODE** button first.
- Press the search buttons  $\blacktriangle$ ,  $\blacktriangledown$  to select
  - **CALIBRATION** mode
- Press **OK**.
- The display shows **PASS 1** - first password number prompt.
- Set with search buttons  $\blacktriangle$ ,  $\blacktriangledown$  the first password number. Press **OK** to confirm it.
- A prompt for the second password number appears on the display - **PASS2**. Enter the second password number in the same way as the first one. Enter the third password number when a **PASS 3** prompt appears on the display.
- Notes: In case an incorrect password is entered, a message WRONG PASSWORD appears on the display. Then make a fresh start.**
- If the right password numbers are entered the display will show **FAT COW MILK**.
- Press the search buttons  $\blacktriangle$ ,  $\blacktriangledown$  to select the milk parameter to be calibrated:
  - **FAT COW MILK** - cow milk fat calibration;
  - **FAT SHEEP (BUFFALO, GOAT) MILK** - sheep (buffalo, goat) milk fat calibration;
  - **SNF COW MILK** - cow milk Solids non fat (SNF) calibration;
  - **SNF SHEEP (BUFFALO, GOAT) MILK** - sheep (buffalo, goat) Solids non fat (SNF) calibration.
  - **DEN COW MILK** - cow milk density calibration;
  - **DEN SHEEP (BUFFALO, GOAT) MILK** - sheep (buffalo, goat) density calibration;
  - **PROT COW MILK** - cow milk protein calibration;
  - **PROT SHEEP (BUFFALO, GOAT) MILK** - sheep (buffalo, goat) protein calibration.
  - **FP COW MILK** - cow milk freezing point calibration;
  - **FP SHEEP (BUFFALO, GOAT) MILK** - sheep (buffalo, goat) freezing

point calibration.

Select desired calibration parameter. Confirm it with **OK**.

- The inscription **VALUE** appears on the display.

- Press the buttons  $\uparrow$ ,  $\downarrow$  to set the calibration value required.

- The number can range from -2,54 to +2,54 (at interval of 0,02) for fat, protein and Solids non fat (SNF), from -12,7 to +12,7 (at interval of 0,10) for density and from -0,250 to + 0,250 (at interval of 0,002) for freezing point calibration. This number defines the calibration value. When desired value is set on the display press the **OK** button to save it.

- The display shows **TOTAL** and a number equal to the total calibration value for this milk parameter. The total calibration value is an algebraic sum of all calibration values added for this milk parameter.

- The calibration value saving is completed.

**Note: The total calibration value can range from -2,54% to +2,54% (for fat, protein and Solids non fat (SNF)), from -12,7°A to +12,7°A for density) and from -0,250 to + 0,250 for freezing point . In case a bigger number has been achieved during a calibration procedure a CAL OUT OF RANGE error message appears on the display.**

---

## SOME REASONS FOR EKOMILK ACCURACY AND REPEATABILITY DEVIATIONS

---

Below are described some reasons that may worsen the accuracy and repeatability of the EKOMILK Analyzers measurement results. Some information about the way to avoid these problems is provided as well.

**1. Aired milk** - this is milk with a lot of tiny air bubbles inside. These air bubbles are very small and it takes long time - from one to more than 10 or even 20 hours for these air bubbles to get out of the milk. This time depends on the milk parameters and mainly on the milk Fat contents - the higher milk fat contents the longer time is required for the air bubbles to get out of the milk. The ultrasonic method is not suitable for aired milk testing since the measurement results are with significant deviations from the real values and even in some particular cases the measurement can not be completed successfully.

The milk becomes aired usually during the milk processing - milking, homogenization, UHT etc. but it can be aired even when the milk sample is mixed if this is made by hard continuous shaking. This is why the sample should be mixed smoothly and carefully.

There are two known methods that allow aired milk to be quickly recovered. The first method requires the measuring mug with the milk sample to be processed for about 10 to 15 seconds in an ultrasonic cleaning machine. The cavitation phenomena of the powerful ultrasonic field removes the air bubbles almost instantaneously. The other method requires the milk sample to be put under pressure - several kg/cm<sup>2</sup> for about 10 to 20 seconds. Unfortunately as you may note both methods require additional equipment but as long as there are small inexpensive ultrasonic cleaning machines commercially available it is easier to apply the first method if necessary.

**2. Milk acidity** - The ultrasonic method requires the milk sample to be warmed during the measurement process. In case the milk has an increased acidity a milk coagulation may occur when the milk is warmed and the measurement results will have significant deviations from the real values and even in some particular cases the measurement can be not completed successfully. For your reference the milk acidity of the milk sample must be less than 25°T for cow, buffalo and goat milk and less than 28°T for sheep milk. Another issue is the milk coagulation contaminates the ultrasonic sensor and it must be properly cleaned then.

**3. Separated milk fat** - If you try to test cool (refrigerated) milk which has some milk fat/cream separated you will probably get wrong result especially for the milk fat contents. In this case you need to warm up the milk up to 40 - 42°C first, mix the milk in order to solve the separated fat, cool it down to 20 - 25°C and then you can test it with the EKOMILK Analyzer.

**4. Contaminated milk** - Any solid particle with a size above 0.5mm may cause measurement result deviations. This is why we recommend the milk sample to be filtered before it is tested if there is a doubt the milk is contaminated.

**5. Milk preservatives** - The milk preservatives change the measurement results. Usually the result deviation is not big but it'll depend on the particular preservative used.

**6. Adulterated milk** - The measurement results may significantly differ from the real milk parameters if the milk contains some additives - salt, sugar, urea etc.

**7. Contaminated sensor** - During the normal work of the EKOMILK analyzers some solid deposits are laid on the ultrasonic sensor walls. In case the analyzer is not regularly and properly cleaned these deposits are gradually accumulated and the measurement results begin to differ from the real milk parameters. This is why it is very important the milk analyzers to be always properly cleaned in accordance with their cleaning procedure.

**8. Power supply** - The power supply can also cause problems with the measurement results accuracy and repeatability. Generally this can happen if the power supply voltage is out of the specified range (220/110V +10 -15%) or if the power supply line is too noisy - especially if there is a powerful equipment working nearby and connected to the same power supply line.

## ERROR LIST

ERROR MESSAGE	PROBLEM & CAUSES & REMEDY
<b>MOTOR ERROR</b>	Motor damage. Contact your dealer to make the repairs. Meanwhile you may use the Analyzer in <b>MOTOR OFF</b> mode. (see “MILK ANALYSIS WITHOUT SUCTION PUMP”).
<b>EMPTY CAMERA ERROR 02</b>	<ol style="list-style-type: none"> <li>1. Air bubble in the measuring sensor. Fill again the mug with the same milk sample and make a fresh start.</li> <li>2. The suction system has some air leak. Check the Plastic Plug position and fix it if it is necessary. If the Plastic Plug stays loose in the Syringe, replace the Syringe. If the problem still exists, test the analyzer in <b>MOTOR OFF</b> mode. If no problem occurs in this mode, the Pump is out of order. Contact your dealer to make the repairs. Meanwhile you may use the Analyzer in <b>MOTOR OFF</b> mode. (see “MILK ANALYSIS WITHOUT SUCTION PUMP”).</li> <li>3. If the problem still exists, the measuring system is contaminated or damaged. Contact your dealer to make the repairs.</li> </ol>
<b>HOT SAMPLE</b>	<ol style="list-style-type: none"> <li>1. The milk sample temperature is +above the spec. Fill again the mug with milk sample with proper temperature and make a fresh start.</li> <li>2. If the problem still exists, the measuring system is contaminated or damaged. Contact your dealer to make the repairs.</li> </ol>
<b>WRONG PASSWORD</b>	Incorrect Password. Reenter the password
<b>CAL OUT OF RANGE</b>	Attempt to enter a calibration value out of allowed limits. See “ <b>CALIBRATION VALUES SAVING</b> ”. If a calibration value out of allowed limits is really needed, the measuring system is out of order. Contact your dealer to make the repairs.
<b>FW MISHMACH</b>	The Main processor is damaged . Contact your dealer for repairs.
<b>ERROR 03</b>	PROM is damaged or erased. Contact your dealer for repairs.
<b>ERROR 07</b>	There is no communication between the main processor and PROM - PROM is burnt or not properly put in the socket. Contact your dealer to make the repairs.
<b>ERROR 09</b>	<ol style="list-style-type: none"> <li>1. The power supply voltage is below spec. Use proper power supply voltage value.</li> <li>2. The measuring system is out of order. Contact your dealer to make the repairs.</li> </ol>
<b>ERROR15</b>	PROM reading/writing problem. Switch off/on the analyzer and try again. If the problem still exists, the PROM is damaged. Contact your dealer to make the repairs.

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