

EKOMILK[®] AMP

Complete system for accurate mastitis and milk quality control

Hybrid Milk Analyzer

EKOMILK HORIZON UNLIMITED (EHU)



EKOMILK[®]
Horizon

Somatic Cell Count (SCC) and milk parameters measurement

- Automated and reliable;
- Fast and cost-effective;
- High accuracy;
- Integrated Wi-Fi module for remote monitoring



ANIMAL MONITORING PLATFORM

AMP

Cloud-based bulk tank, dairy herd & cow health management and instrument maintenance
www.animalmonitoring.com



SAMPLE FEEDER

Peripheral device for automatic milk sampling and feeding



- Flexible process automation;
- Autonomous processing of 49 milk samples;
- Processing time around 1 hour;
- Auto cleaning function

approx.
€0.03
per test



Certificate delivered by the European Commission, as the institution managing Horizon 2020, the EU Framework Programme for Research and Innovation 2014-2020

The project proposal 756888, MastCloud

System Applications:

- ◆ dairy farms with a conventional milking parlor ◆
- ◆ dairy farms with robotic milking machines ◆
- ◆ large laboratories testing over 100 000 samples / year ◆
- ◆ small labs ◆ vet cabinets ◆ milk processors ◆
- ◆ mobile labs for milk collection ◆

EKOMILK[®] AMP

SOMATIC CELL COUNT & MILK PARAMETERS

Milk parameters	Cow / Tank	Measuring Range and Accuracy
Somatic Cell Count	Cow / Tank	80 000 ÷ 15 000 000 cells/ml ± 8%
Fat	Cow / Tank	0,5% ÷ 12% ± 0,04%
Protein	Cow / Tank	1,0% ÷ 9,0% ± 0,14%
Solids not fat (SNF)	Cow / Tank	6% ÷ 12% ± 0,2%
Temperature	Cow / Tank	0°C ÷ 50°C ± 0,1%
Conductivity	Cow / Tank	1,0 mS/cm ÷ 10 mS/cm ± 0,2 mS/cm
Density	Tank only	1,0200 ÷ 1,0400 g/cm ³ ± 0,0005 g/cm ³
Freezing point	Tank only	-1,000°C ÷ 0,000°C ± 0,015°C
Added water	Tank only	0,0% ÷ 60% ± 5%
Lactose	Cow / Tank	0,5% ÷ 7% ± 0,2%

- ◆ The accuracy of SCC, Protein, Fat and SNF are achieved after AMP calibration finetuning and valid for ICAR-compliant milk samplers;
- ◆ Lactose indicator is available only via AMP on demand;
- ◆ All specifications are preliminary and can be changed without prior notice

SYSTEM CONFIGURATIONS



EKOMILK[®] AMP

Mini-lab for fully-automated, on-farm, fast and reliable SCC and milk parameter analysis at €0.03/test in around a minute per measurement

Animal Monitoring Platform - Result Transfer, Workflow Integration and Measurement Supervision

AMP Standard package:

- Admin User Account with up to 1000 cows per farm;
- Automatic Data Transfer via integrated solid Wifi and USB stick workflow support;
- Powerful online calibration engine which allows to obtain unparalleled accuracy over an extended range of milk parameters;
- Automatic dashboards on bulk tank, on individual cow and on herd level;
- Real-time raw results viewing; mastitis visualization and monitoring;
- Efficient sample identification processes;
- Instrument maintenance program;
- Comprehensive notification platform via SMS, email and apps - receive measurement results and instrument warnings on any device anywhere

Optional:

- Up to 20 000 cows per farm and unlimited farms;
- Customized calibration services incl. automatic seasonal calibration update;
- Extra animal health indicators - lactose, acidosis and ketosis.
- Custom branding admin interface & 3rd party modules & APIs

1-year standard AMP admin account is included with every system

EKOMILK[®] Horizon

ESSENTIAL (EHE) +AMP

Automated milk analyzer for fast and effective control of SCC in dairy farms and processing plants

EKOMILK[®] Horizon

UNLIMITED (EHU) +AMP

Automated hybrid milk analyzer for SCC and milk parameters measurement

EKOMILK[®] Horizon

ROBOT VERSION (EHR) +AMP

Direct simple connection to milking robots; Real-time analysis allowing immediate action for sick animals

EKOMILK[®] AMPI SAMPLE FEEDER (SF) +AMP

Peripheral device for automatic milk sampling and feeding to the analyzer for fast milk parameters' analysis

ISO 9001
CE

Buy online at www.ekomilkhorizon.com



BULTEH 2000 Ltd.
18, Promishlena Str.,
Industrial zone 19,
6000, Stara Zagora
BULGARIA

tel./fax: (+359) 42 626019
tel.: (+359) 42 220883
e-mail: bulte@bulteh.com; info@eonbg.com
web: www.bulteh.com
e-commerce: www.ekomilkhorizon.com

AMP Innovation Ltd.
e-mail: office@animalmonitoring.com
web: www.ekomilk.eu
web: www.animalmonitoring.com